


MIRA Ophthalmic Cryo

An ISO 13485 Company



*The Ophthalmic Cryo
that provides instantaneous
freezing and defrosting
with tip-freeze-only
probes*

 *Approved*



Optimized Design Means Less Tissue Damage

Slow freezing and refreezing produce large ice crystals that destroy cell membranes and cause excessive tissue damage. Under such conditions, inadvertent freezing of extraocular tissues, including the lids, may cause unnecessary post-operative patient discomfort, depigmentation, and alopecia. Both the rate of freezing and thawing are important. In fact, slow thawing causes even more tissue damage than slow freezing. Both increase the size of the ice crystals in the tissues. Fast freezing and thawing, however, produce delicate intracellular crystal formation. This is less destructive and is useful in treating retinal detachment. Thus a very fast freeze and defrost cycle is preferable to a slower one.¹ For this fundamentally significant reason, MIRA offers the fastest freeze and defrost system available. The MIRA Ophthalmic Cryo is the only cryo that provides instantaneous freezing and defrosting with all tip-freeze-only probes. This is because the MIRA Cryo is the only unit designed and optimized specifically for ophthalmic use. Our instant freeze and defrost feature also allows you to operate more quickly, thus improving your efficiency and productivity.

MIRA's exclusive design, based on the Joule-Thompson principle, utilizes a high pressure system in both the freeze and defrost modes. This system guarantees instantaneous freezing and defrosting. When the footswitch is depressed, high pressure gas is forced along the narrow periphery, through a narrow orifice, and into a larger volume chamber in the tip. This sudden expansion of the gas causes an instantaneous heat loss, thereby abruptly freezing the probe tip. The console temperature selector automatically regulates the gas pressure. This allows a preset frozen tip temperature to be maintained in direct contact with the tissue to be treated.

Our unique rapid defrost occurs when the footswitch is released and the high pressure gas is forced in the reverse direction. This sudden compression of the gas causes an instantaneous gain in heat, instantly warming the tip.



- **Better Visibility and Maneuverability**
- **Probes For All Applications**
- **No Manual Pressure Adjustments Required**
- **Safe, Quiet, Reliable Operation**
- **Fast Installation and Disassembly**

References

1. Pruett, R.C. Ophthalmic Cryo, *The FORUM*, MIRA, Inc., Waltham, MA, 1988, Vol. 1, No. 3, pp. 1, 3.

An ISO 13485 Company

Makes the Surgeon's Task Easier

MIRA has designed its Cryo with the surgeon's needs in mind. Our tubing with *enhanced flexibility and natural downward positioning of the probe tip* are intended to make the vitreoretinal specialist's work easier. Our probe tubing is a full nine feet long in order to provide the surgeon with a *greater range of motion*. With *reduced coil memory*, the tubing will not offer resistance when maneuvering the probe or spontaneously recoil when placed upon an operating table. Another benefit of the unique MIRA design is that *the handle remains at a comfortable temperature — even during extended freezing procedures*. This is because the warm incoming gas completely surrounds the central tube of cold exhaust gas, thus insulating the probe handle.

Practical Office System


Our efficient design does not stop at just the cryo console and probes. For added convenience, *both our single and dual carriers can accept either small diameter or large diameter cylinders*. A dual switchover valve on the two cylinder carrier allows continued operation of the cryo while an empty cylinder is removed for refill.

Durable Construction

MIRA has taken extra steps to ensure the unit's expected long life. All probe handles are made of stainless steel. Both probe tips and connector jacks come equipped with protective covers. Because of the corrosive nature of CO₂ and N₂O gases, our gas cylinder pressure lines are made of *copper tubing to reduce the risk of damage* to these lines. *Replaceable microfiber filters* can be fitted above the gas cylinder connections. These filters are highly recommended, especially for high humidity climates and countries where gas contamination is a problem. The cryo filter removes condensed water, oil droplets, other liquids, and solid particles from the gas. Using a filter will help *extend the life of your probes* by preventing these contaminants from clogging the gas flow passages.

One Year Warranty

MIRA manufactures only the highest quality products, and our Cryo is no exception. The MIRA Cryo comes with a one year warranty on both materials and workmanship. If a repair is ever needed, we will respond with prompt, efficient service.





CR4010 Curved Retinal Probe
2.8 mm diameter x 17.3 mm length



CR4022 Straight Cataract Probe
2.1 mm diameter x 27.0 mm length



CR4023 Curved Cataract Probe
2.1 mm diameter x 19 mm length



CR4025 Mini Curved Cataract Probe
1.5 mm diameter x 19.0 mm length



CR4030 Curved Glaucoma Probe
3.4 mm diameter x 19.0 mm length



CR4040 Vitreous Probe
1.5 mm diameter x 27.0 mm length



CR4045 Mini Vitreous Probe
1.0 mm diameter x 31.6 mm length



CR4075 Hammerhead Probe
2.8 mm diameter x 26.5 mm length
6 mm tip width



CR4080 Baby Hammerhead Probe
1.5 mm diameter x 22.5 mm length
4 mm tip width



CR4085 Gaynon ROP Probe
2.0 mm diameter x 19.0 mm length
5 mm freezing length



CR4090 Retinoblastoma Probe
1.5 mm diameter x 35.0 length
10 mm radius



Cryo Systems

- CR4003 Ophthalmic Cryo O.R. System Includes:**
- CR4000 Ophthalmic cryo console with footswitch
 - CR4010 (2) 2.5 mm curved retinal probe
 - CR4150 Dual cylinder Cryo carrier
 - CR4160 25-foot exhaust hose
 - CR5350 Sintered Metal Filters (pack of 3)
 - CR43xx Gas connector/filter (specify type of gas)
 - CR40xx Cryo Instructional Video (Specify PAL or SECAM)
- CR4001 Ophthalmic Cryo System Includes:**
- CR4000 Ophthalmic cryo console with footswitch
 - CR4010 2.5 mm curved retinal probe
 - CR4160 25-foot exhaust hose
 - CR4260 60 inch pressure hose connection
 - CR4350 12 Filters and 6 "O" Rings
 - CR43xx Gas connector/filter (specify type of gas)
 - CR40xx Cryo Instructional Video (Specify PAL or SECAM)

Cryo Specifications

Console Dimensions:	Width 5 inches (13.3 cm) Height 6.5 inches (16.5 cm) Depth 11 inches (27.9 cm)
Weight:	16 pounds (7.3 kg)
Front Panel Gauge:	Indicates incoming cylinder gas pressure
Front Panel 3 Position Temperature Selector:	-25°C, -55°C, -85°C, tolerance $\pm 5^\circ\text{C}$
Front Panel On/Off Switch:	Turns console on/off (controls gas lines)
Footswitch:	Controls freezing operation (depress to freeze, release to defrost)
Power Source:	Operates on CO ₂ or N ₂ O gas (no electrical power required)

Cryo Carriers

- CR4150** Dual cylinder cryo carrier with three 30 inch hoses & one switchover valve (holds either 2 C or 2 E cylinders)

Cryo Filters

- CR4400 Microfiber Filter Retrofit
- CR4350 12 Filters and 6 "O" Rings
- CR5350 Sintered metal Filters (pack of 3)
- CR4148 Filter retro-fit kit for CR4150
- CR5400 Sintered metal filter

Gas Connectors

- CR4210 Handtight CGA N₂O for C cylinder
- CR4220 Handtight CGA CO₂ for C cylinder
- CR4250 Universal T-yoke N₂O or CO₂ for E cylinder
(available for export only)
- CR4280 T-Yoke N₂O for E cylinder
- CR4290 T-Yoke CO₂ for E cylinder
- CR4310 Handtight CGA N₂O for C cylinder with filter
- CR4320 Handtight CGA CO₂ for C cylinder with filter
- CR4380 T-Yoke N₂O for E cylinder with filter
- CR4390 T-Yoke CO₂ for E cylinder with filter
- CR5310 Handtight N₂O for "C" cylinder with sintered metal filter
- CR5320 Handtight CO₂ for "C" cylinder with sintered metal filter

Cryo Accessories

- CR4160 25 foot exhaust hose
- CR4260 60 inch pressure hose connection
- CR4360 30 inch pressure hose
- CR4270 Protective probe tip cap